	•	26, 2000 v. None				
Batavia, IL 60510 CMS ME1/2 A ELECTRIC	elerator Laboratory NODE PANEL AL TESTING VELER					
Reference	Drawing(s)					
Endcap Muon Chambe	er ME1/2 Final Assembly					
<u> -</u>	E-368120					
Endcap Muon Chamber Anode Panel Assembly						
	E-368121					
Budget Code:	Project Code:					
Released by:	Date:					
Prepared by: M. Hubbard, B. Jensen, L. Lee						
Title	Signature D	ate				
TD / E&F Process Engineering	Bob Jensen/Designee					
TD / E&F CMS Assembly	Glenn Smith/Designee					
TD /E&F Technological Physicist	Oleg Prokofiev/Designee					
TD / E&F CMS Project Manager	Giorgio Apollinari/Designee					

Panel Serial No.

Revision Page

			<u>Revision Page</u>		
Revision	Step No.		Revision Description	TRR No.	Date
None	N/A	Initial Release		N/A	04/26/00
		IHEP			
		IHEP			

Ensure appropriate memos and specific instructions are placed with the traveler before issuing the sub traveler binder to production.


1.0	Gener	al Notes
	1.1	White (Lint Free) Gloves (Fermi stock 2250-1800) or Nitrile Gloves (Fermi stock 2250-2040) shall be worn by all personnel when handling all product parts after the parts have been prepared/cleaned.
	1.2	All steps that require a sign-off shall include the Technician/Inspectors first initial and full last name.
	1.3	No erasures or white out will be permitted to any documentation. All incorrectly entered data shall be corrected by placing a single line through the error, initial and date the error before adding the correct data
	1.4	All Discrepancy Reports issued shall be recorded in the left margin next to the applicable step.
	1.5	All personnel performing steps in this traveler must have documented training for this traveler and associated operating procedures.
	1.6	Personnel shall perform all tasks in accordance with current applicable ES&H guidelines and those specified within the step.
	1.7	Cover the product/assembly with Mylar when not being serviced or assembled.
	1.8	Never hand/pass anything over a panel as dropped items may damage the panel.
2.0	Parts 1	Kit List
	2.1	Attach the completed Parts Kit for this production operation to this traveler. Ensure that the serial number on the Parts Kit matches the serial number of this traveler. Verify that the Parts Kit received is complete.
		Process Engineering/Designee Date
		<u> </u>

Panel Serial No._____

3.0	Panel Pr	eparation		Rev. None
	3.1	Acquire the Anode Panel (ME-368121) as per to of this traveler. Put the Anode panel on transport		Completed
	3.2	Rotate the panel until in the vertical position. Introgen hand airgun (MX-????) to remove an on the panel		
	Note(s):	Extreme care must be used while using the Ior to prevent damage to the Anode Panel Wires.		
	3.3	Safety Glasses must be worn while using the I Inspect the Anode panel (both sides) to be sur materials removed from the panel surface.		
		Technician(s)	Date	

4.0 <u>Cathode Strip Resistance tests</u>

4.1 Using a Multimeter, and a Toggle Switch Box, check the continuity in resistance of the cathode strip connectors. Starting with the connector left of the serial number, test each connector and if it passes, check it off in the chart below. If it fails, write the




Note(s):

All measurements must be within the range of 0.9 – 1.1 Meg Ohm.

	Pass	Fail	
Connector #1			
Connector #2			TT)
Connector #3			
Connector #4			
Connector #5			

Remarks:			

Note(s):		
	After measurements are complete	d inform supervisor of any failure
	If all strips pass, panel is acceptable	le to continue.
	_	_



CMS ME1/2 Anode Panel Electrical Testing

5.0 Anode Wire Group Capacitance Measurements

Using a Capacitance Measuring Unit, measure the anode wire group capacitance from the protection boards. Measuring will be start from the narrow side of panel.

Note(s): After measurements are completed inform supervisor of any discrepancy with reference data table with the capacitance measurements. Be sure that cathode connectors on the box are terminated to ground. Protection boards 7 8


		Protection Board								Replaced
		1	2	3	4	5	6	7	8	
	1									
C	2									
H	3									
A	4					TI		<u> </u>	T 1	
N	5					_				
N	6									
E	7							<u> </u>		
L	8									
	9									
N	10									
U	11									
M	12									
В	13									
Е	14					4				
R	15		r-7			-7		r.	- 1	
	16									
Rai LOW	nge: ?HIG									
	Remark	xs:	<u> </u>	1	1	1		, r	1	
										-

CMS ME1/2 Anode Panel Electrical Testing

Technician(s) Date

6.0 <u>Strip to Ground Capacitance Measurement</u>

6.1 Using a switch box, cable and LCR meter, measure the Capacitance from Strip to Ground.



		Cathode Connector							
		1	2	3	4	5			
	1								
	2								
C	3								
H	4								
A N	5		D.						
N	6								
E	7		7						
L	8								
	9								
	10								
N	11								
U M	12								
B	13		1			1			
E	14								
R	15	T			T				
	16								
Ra LOV	nge: V?HI								
	Remar	ks:							
			7						
						\dashv			

CMS ME1/2 Anode Panel Electrical Testing

Technician(s) Date

CMS ME1/2 Anode Panel Electrical Testing

Completed


7.0	High	Voltage	Tests

7.1 Transport the Anode panel to the High Voltage Test station.

7.2 Install the Anode Panel into the High Voltage Test Box. In accordance with below diagram, connect test cables to the following:

Protection Boards (Anode Cables #1 - #8)
Cathode connectors (Cathode Cables #1 - #5)

High Voltage banana plugs (High Voltage Cables #1 - #3)



Technician(s) Date

r








7

7.3 Close the box, tighten the seal bolts and open valve with dry air. Purge dry air through the box for approximately 10 - 15 minutes until the humidity will drop to 30% or less. Record the start humidity and the ending humidity.



7.4 Switch on the High Voltage Power Supply. Slowly raise the High Voltage up to 1.0 kV and allow the value to stabilize for 2-3 minutes. After the value is stable, record it on the chart below. Continue to increase the voltage to the next HV point, allowing the value to stabilize for 2-3 minutes each time, and record the data for the values onto the chart up to and including 4.5 kV.

High V	Strip Side	Test #1	Ļ		Non-Strip Side		4	
HV (kV)	I (mA)	Start Time	Stop Time	Humidity %	I (mA)	Start Time	Stop Time	Humidity %
1.0								
2.0								
3.0								
3.5		r	1	r	1			
4.0								
4.1								
4.2					4			
4.3			r	Γ	T			
4.4								
4.5								
Remarks	:							

Note(s):

If a problem arises, lower the kV to ZERO.

All measured current data will be lower are completed inform supervisor of any		
Technician(s)	Date	

CMS ME1/2 Anode Panel Electrical Testing

High Voltage: Test #2

7.5 Perform High Voltage: Test #2 only if discrepancies occurred in Step 6.4 and repairs were required.

Non-Strip

	Side			Side				
HV	I (mA)	Start Stop	Humidity %	I (mA)	Start	Stop	Humidity %	
(kV)		Time Time		- 1	Time	Time		
1.0								
2.0								
3.0								
3.5								
4.0						r		
4.1								
4.2								
4.3				4				
4.4				T				
4.5								
Remarks	3:							
Note(s):	Note(s):	a problem arises, lowe	er the kV to ZERO.					
All measured current data will be lower than 1 μA . After all measurements are completed inform supervisor of any discrepancy with reference data table.								
	Τε	echnician(s)			Date			
CMS MI	E1/2 Anode	Panel Electrical Testin	g		Panel Ser	ial No		

7.6 Keep plane under $HV=4.5kV \pm .005kV$ for about 3-5 minutes. Measure current from each HV segment at $4.5 \, kV \pm .005kV$. Record data into below table.

Segment Number							
	All 1 Segments		2	3	4	5	Humidity
Plane #1	μΑ	μА	μA	μА	μΑ	μΑ	%
Strip Side						T	
Plane #2	μΑ	μΑ	μΑ	μΑ	μΑ	μΑ	%
Non-Strip Side							
Note(s):	All measured current of						
	with reference data tab				Date		
	(-)						Completed
7.7	If the current in any seg air or Ethyl Alcohol (Fe 2500). Then repeat test	rmi Stk. No.	. #1920-0600) and	d a low-lint w	vipe (Fermi Stk		
		Ses	gment Numb	er		T]	
	All Segments		2	3	4	5	Humidity
Plane #1 Strip Side	μA	μΑ	μА	μΑ	μΑ	μΑ	%
Plane #2 Non-Strip Side	μA	μΑ	μА	μА	μA	μА	%
Note(s): All measured current data from segment will be lower than 0.1 $\mu A.$							
7.8	Decrease the HV slowly Close air gas manifold.	y to 0 Volts	and switch off H	ligh Voltage	Power Supply		
7.9	Loosen the bolts and re	emove anod	e panel.	7			
7.10	Install panel on the par	nel cart and	transport panel t	to the storage	e area.		
	Technician(s)				Date		

Panel Serial No.

CMS ME1/2 Anode Panel Electrical Testing

Production Complete 8.0

	8.1	Process Engineering verify that the CM accurate and complete. This shall inclu completed and signed off. Ensure that Forms, Deviation Index and disposition before being approved.	ade a review of all st all Discrepancy Rep	teps to ensure theorts, Nonconfo	hat all operations have rmance Reports, Repair	been Rework
		Comments:				
		Process Engineering/Designee		Date		
.0	Attach	the Process Engineering O.K. to Proceed	Tag to the device th	nis production (operation performed.	
		Process Engineering/Designee	TT	Date		
0.0	Procee	d to the next major assembly operation.				
		Process Engineering/Designee		Date		
		7 1	r 1		1	
MS I	ME1/2 An	node Panel Electrical Testing		Panel Se	rial No	

Page 14 of 14